

To select the correct size, measure your leg from the tip of the toe to the top of the AFO / KAFO. Select the sock which is closest to that length.

There is no heel pocket, so the shoe size doesn't matter. The sock is like a long stretchy tube that will mould to the contours of the foot and leg. Measurements are the unstretched socks.

SmartKnit AFO Inte	rface Seamless Socks	
Size	Width	Length
Infant Small AFO	3.5cm	18cm
Infant AFO	5cm	23cm
Child Small AFO	5.5cm	28cm
Child AFO	6.5cm	36cm
Child Extra Long AFO	7.5cm	51cm
Adult Small AFO	7.5cm	46cm
Adult Regular AFO	9cm	51cm

Patterned AFO Core-Spun Interface Socks				
Size	Width	Length		
Child Small	5.5cm	28cm		
Child Regular	6.5cm	36cm		
Child Extra-Long	7cm	51cm		
Adult Small	7.5cm	46cm		
Adult Regular	9cm	51cm		

Website:	Company No.:	Contact Details:	Address:
www.sensorysmart.co.uk	09477873	help@sensorysmart.co.uk	Concorde House, Union Drive, Sutton Coldfield,
			West Midlands, B73 5TE, UK



SmartKnit Big Toe AFO Seamless Socks			
Size	Width	Length	
Infant	5cm	23cm	
Child Small	5.5cm	28cm	
Child Regular	6.5cm	36cm	
Adult Small	8cm	36cm	
Adult Regular	9cm	51cm	

SmartKnit KAFO Interface Seamless Socks

Size	Width	Length
Child XX-Long KAFO	6.5cm	46cm
Adult Narrow KAFO	10 x 6.5cm	61cm
Adult Regular KAFO	15 x 9cm	71cm

Torso Interface® — Turtleneck with Sleeves

Size	Width	Length	Fits Circumference
X-Small	4" (10cm)	12" (30cm)	10"-12" (25-30cm)
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)

Website:	Company No.:	Contact Details:	Address:
www.sensorysmart.co.uk	09477873	help@sensorysmart.co.uk	Concorde House, Union Drive, Sutton Coldfield,
			West Midlands, B73 5TE, UK



Torso Interface® — Without Straps

Size	Width	Length	Fits Circumference
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)
Medium Long	8" (20cm)	20" (50cm)	18"-24" (45-60cm)
Large	10" (25cm)	20" (50cm)	22"-30" (55-75cm)
Large Long	10" (25cm)	24" (60cm)	22"-30" (55-75cm)
X-Large	12" (30cm)	24" (60cm)	26"-38" (66-95cm)
KS	14" (35cm)	24" (60cm)	34"-48" (76-122cm)

Torso Interface® — Crew-Neck

Size	Width	Length	Fits Circumference
X-Small	4" (10cm)	12" (30cm)	10"-12" (25-30cm)
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)
Medium Long	8" (20cm)	22" (56cm)	18"-24" (45-60cm)
Large	10" (25cm)	20" (50cm)	22"-30" (55-75cm)
Large Long	10" (25cm)	24" (60cm)	22"-30" (55-75cm)
X-Large	12" (30cm)	24" (60cm)	26"-38" (66-95cm)
KS	14" (35cm)	24" (60cm)	30"-48" (76-122cm)
кх	14" (35cm)	28" (71cm)	30"-48" (76-122cm)

Website:	Company No.:	Contact Details:	Address:
www.sensorysmart.co.uk	09477873	help@sensorysmart.co.uk	Concorde House, Union Drive, Sutton Coldfield,
			West Midlands, B73 5TE, UK



KnitRite Seamless AFO, KAFO & Orthotic Size Chart

Torso Interface® — Crew-Neck: Axilla Flaps

Size	Width	Length	Fits Circumference
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)
Medium Long	8" (20cm)	22" (56cm)	18"-24" (45-60cm)
Large	10" (25cm)	20" (50cm)	22"-30" (55-75cm)
Large Long	10" (25cm)	24" (60cm)	22"-30" (55-75cm)
X-Large	12" (30cm)	24" (60cm)	26"-38" (66-95cm)

Torso Interface[®] — Crew-Neck: With Sleeves

Size	Width	Length	Fits Circumference
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)
Medium Long	8" (20cm)	22" (56cm)	18"-24" (45-60cm)
Large	10" (25cm)	20" (50cm)	22"-30" (55-75cm)
Large Long	10" (25cm)	24" (60cm)	22"-30" (55-75cm)
X-Large	12" (30cm)	24" (60cm)	26"-38" (66-95cm)
XX-Large	12" (30cm)	28" (71cm)	26"-38" (66-95cm)

Website:	Company No.:	Contact Details:	Address:
www.sensorysmart.co.uk	09477873	help@sensorysmart.co.uk	Concorde House, Union Drive, Sutton Coldfield,
			West Midlands B73 5TE LIK



KnitRite Seamless AFO, KAFO & Orthotic Size Chart

Size	Width	Length	Fits Circumference
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)
Medium Long	8" (20cm)	22" (56cm)	18"-24" (45-60cm)
Large	10" (25cm)	20" (50cm)	22"-30" (55-75cm)
Large Long	10" (25cm)	24" (60cm)	22"-30" (55-75cm)

Torso Interface® — V-Neck: Left Flap

Size	Width	Length	Fits Circumference
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)
Large	10" (25cm)	20" (50cm)	22"-30" (55-75cm)
Large Long	10" (25cm)	24" (60cm)	22"-30" (55-75cm)
X-Large	12" (30cm)	24" (60cm)	26"-38" (66-95cm)

Torso Interface® — V-Neck: No Flaps

Size	Width	Length	Fits Circumference
Small	6" (15cm)	14" (35cm)	14"-18" (35-45cm)
Medium	8" (20cm)	16" (40cm)	18"-24" (45-60cm)
Medium Long	8" (20cm)	20" (50cm)	18"-24" (45-60cm)
Large	10" (25cm)	20" (50cm)	22"-30" (55-75cm)
Large Long	10" (25cm)	24" (60cm)	22"-30" (55-75cm)
X-Large	12" (30cm)	24" (60cm)	26"-38" (66-95cm)
KS	14"(35cm)	24" (60cm)	30"-48" (76-122cm)

Website:	Company No.:	Contact Details:	Address:
www.sensorysmart.co.uk	09477873	help@sensorysmart.co.uk	Concorde House, Union Drive, Sutton Coldfield, West Midlands, B73 5TE, UK