

## Spio Size Chart


### SPIO TLSO System (Vest)

Choose size based on chest measurement.  
 Take measurement with patient lying down.

**Available Sizes:**

V 41 V 44 V 47 V 51 V 54  
 V 60 V 67 V 70 V 73 V 77

Chest  
 circumference in  
 centimetres at  
 nipple line




### SPIO Upper Body Orthosis (Shirt)

Choose size based on chest measurement.  
 Take measurement with patient lying down.  
 Short and long sleeve available.

**Available Sizes:**

S 41 S 44 S 47 S 51 S 54 S 57  
 S 60 S 64 S 67 S 70 S 73 S 77

Chest  
 circumference in  
 centimetres at  
 nipple line




### SPIO Lower Body Orthosis (Pant)

Choose size based on waist measurement.  
 Take measurement with patient standing up.  
 Short and long pant available.

**Available Sizes:**

P 41 P 43 P 45 P 47 P 49 P 51  
 P 53 P 55 P 57 P 59 P 61 P 63  
 P 65 P 67 P 69 P 71 P 73

Waist  
 circumference in  
 centimetres at  
 navel



### SPIO Wrist Hand Orthosis (Glove)

Choose size based on palm measurement.

**Available Sizes:**

G 4 G 4.5 G 5 G 5.5 G 6 G 6.5  
 G 7 G 7.5 G 8

Measure across  
 palm in centimetres

